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22.2 Funkmodul (Bluetooth)

PAN1322-SPP
ENW89841A3KF

Bluetooth Qualification and Regulatory Certification

ENW89841A3KF is intended to be installed inside end user equipment. ENW89841A3KF is Bluetooth-qualified and also FCC-certified and Industry Canada approved, and conforms to R&TTE (European) requirements and directives with the reference design described in [Figure 9](#).

Manufacturers of mobile, fixed or portable devices incorporating this device are advised to clarify any regulatory questions and to have their complete product tested and approved for compliance (FCC or other when applicable). When using other antennas, a "class II permissive change" is required for FCC approval. The normal procedure is to first provide a technical test report showing that 4 dBi is not exceeded and to continue working with a regulatory test house to finalize the approval for a new antenna implementation.

There are no parts in ENW89841A3KF that can be modified by the user except modifications of the device BD data and loading of SW patches. Any changes or modifications made to this device that are not expressly approved by Panasonic, may void the user's authority to operate the equipment.

9.2 FCC Class B Digital Devices Regulatory Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by 1 or more of the following measures:

- Reorient or relocate the antenna
- Increase the separation between the equipment and receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio or television technician for help

9.3 FCC Wireless Notice

This product emits radio frequency energy, but the radiated output power of this device is far below the FCC radio frequency exposure limits. Nevertheless, the device should be used in such a manner that the potential for human contact with the antenna during normal operation is minimized.

To meet the FCC's RF exposure rules and regulations:

- The system antenna used for this transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
- The system antenna used for this module must not exceed 4 dBi.
- Users and installers must be provided with antenna installation instructions and transmitter operating conditions for satisfying RF exposure compliance, please refer to [Figure 10](#).


Figure 10 Cutout Drawing

Manufacturers of mobile, fixed or portable devices incorporating this module are advised to clarify any regulatory questions and to have their complete product tested and approved for FCC compliance.

9.4 FCC Interference Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference

2. This device must accept any interference received, including interference that may cause undesired operation.

9.5 FCC Identifier

FCC ID: T7VEBMU

9.6 European R&TTE Declaration of Conformity

Hereby, Panasonic Industrial Devices Europe GmbH, declares that the Bluetooth module ENW89841A3KF is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC.

As a result of the conformity assessment procedure described in Annex III of the Directive 1999/5/EC, the end-customer equipment should be labelled as follows:

Figure 11 Equipment Label

PAN1322 In the specified reference design can be used in the following countries:

Austria, Belgium, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Poland, Portugal, Slovakia, Slovenia, Spain, Sweden, The Netherlands, the United Kingdom, Switzerland, and Norway.

Declaration of Conformity (DoC) 1999/5/EC									
We,	<u>Panasonic Industrial Devices Europe GmbH</u> <u>Wireless Connectivity, Power Electronics R&D Center</u> <u>Zeppelinstrasse 19, 21337 Lueneburg, Germany</u>								
declare under our sole responsibility that the product:									
Type of equipment:	<u>Bluetooth Module</u>								
Brand name:	<u>PAN1321 / PAN1311</u> <u>PAN1322 / PAN1312</u>								
Model name:	<u>ENW89811K4CF / ENW89810K5CF</u> <u>ENW89841A3KF / ENW89841C3KF</u>								
to which this declaration relates, is in compliance with all the applicable essential requirements, and other provisions of the European Council Directive:									
1999/5/EC	Radio and Telecommunications Terminal Equipment Directive (R&TTE)								
The conformity assessment procedure used for this declaration is Annex IV of this Directive.									
Product compliance has been demonstrated on the basis of:									
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 2px;">- EN 50371: 2002-11</td> <td rowspan="2" style="padding: 2px;">For article 3.1 (a) : Health and Safety of the User</td> </tr> <tr> <td style="padding: 2px;">- EN 60950-1: 2011-01</td> </tr> <tr> <td style="padding: 2px;">- EN 301 489-1 V1.9.1 (2011-04)</td> <td rowspan="2" style="padding: 2px;">For article 3.1 (b) : Electromagnetic Compatibility</td> </tr> <tr> <td style="padding: 2px;">- EN 301 489-17 V2.1.1 (2009-05)</td> </tr> <tr> <td style="padding: 2px;">- EN 300 328 V1.7.1 (2006-10)</td> <td style="padding: 2px;">For article 3.2 : Effective use of spectrum allocated</td> </tr> </table>	- EN 50371: 2002-11	For article 3.1 (a) : Health and Safety of the User	- EN 60950-1: 2011-01	- EN 301 489-1 V1.9.1 (2011-04)	For article 3.1 (b) : Electromagnetic Compatibility	- EN 301 489-17 V2.1.1 (2009-05)	- EN 300 328 V1.7.1 (2006-10)	For article 3.2 : Effective use of spectrum allocated	
- EN 50371: 2002-11	For article 3.1 (a) : Health and Safety of the User								
- EN 60950-1: 2011-01									
- EN 301 489-1 V1.9.1 (2011-04)	For article 3.1 (b) : Electromagnetic Compatibility								
- EN 301 489-17 V2.1.1 (2009-05)									
- EN 300 328 V1.7.1 (2006-10)	For article 3.2 : Effective use of spectrum allocated								
The technical construction file is kept available at:									
<u>Panasonic Industrial Devices Europe GmbH, Zeppelinstrasse 19, 21337 Lueneburg, Germany</u>									
Issued on:	<u>31st of October 2012</u>								
Signed by the manufacturer:									
(Company name)	<u>Panasonic Industrial Devices Europe GmbH</u>								
(Signature)	<u>Panasonic Industrial Devices Europe GmbH</u> <u>Zeppelinstrasse 19</u> <u>21337 Lueneburg</u> <u>Tel.: +49 (0) 4131 / 899-0</u>								
(Printed name)	<u>Heino Kaehler</u>								
(Title)	<u>Manager Wireless Connectivity</u>								

Figure 12 Declaration of Conformity

9.7 Bluetooth Qualified Design ID

Panasonic has submitted End Product Listing (EPL) for PAN1322, based on Intel eBMU platform, in the Qualified Product List of the Bluetooth SIG. These EPL are referring the Bluetooth qualification of the SPP-AT application running on the eBMU chip under QD ID B021246.

Manufacturers of Bluetooth devices incorporating PAN1322 can reference the same QD ID number.

Bluetooth QD ID: B021246 (PAN1322 SPP BT2.1).

9.8 Industry Canada Certification

PAN1322 complies with the regulatory requirements of Industry Canada (IC), license: IC: 216Q-EBMU

Manufacturers of mobile, fixed or portable devices incorporating this module are advised to clarify any regulatory questions and ensure compliance for SAR and/or RF exposure limits. Users can obtain Canadian Information on RF exposure and compliance from www.ic.gc.ca.

This device has been designed to operate with the built in antenna. It is not allowed to alter the antenna or connecting an external antenna to the module. The built in antenna used for this transmitter must not be collocated or operating in conjunction with any other antenna or transmitter.

9.9 Label Design of the Host Product

It is recommended to include the following information on the host product label:

Contains transmitter Module FCC ID: T7VEBMU / IC: 216QEBMU

9.10 Regulatory Test House

The test house used by Panasonic in the Bluetooth and Regulatory approvals for the module PAN1322:

Eurofins Product Service GmbH
Storkower Str. 38c
D-15526 Reichenwalde b. Berlin
GERMANY
Tel.: +49 33631 888 0
Fax: +49 33631 888 650
www.eurofins.com